ENERGI Advanced Boot Camp Program



Super Set "Partner Up"

Dynamic Warm Up: All 24 participants will complete the following exercises while moving forward approximately 10 yards. Repeat sequence for a total of 2 sets.

- Alternating Heel Tucks
- Walking Lunge
- Inch Worm

- Straight Leg Lift
- Jump Forward (Both Feet)
- Ankle Rotations

Group #1

	RGI Station 3 -Up & Balance, and Squat & Dip
1	Triceps Push-Up (6 each leg)
2	Squat to V-Sit

Group #2

,	Jump Rope Intervals	
	1	Two Foot
	2	"Rocky" Style - Alternating Right to Left Foot

Group #3

ENERGI Station 3 Push-Up & Balance	
1	"T" Push-Up (6 each side)
2	Single Leg "T" (6 each side-tempo)

Group #4

Jui	mp Rope Intervals
1	Two Foot – Front-to-Back
2	"Rocky" Style - Alternating Right to Left Foot

Group #5

	ERGI Station 4 lity & Power Steps
1	Fast Toe Up (15 with right then switch to left foot)
2	Lateral Step-Up (10x right leg then switch 10x left leg)

Group #6

Jump Rope Intervals	
1	Two Foot – Side-to-Side
2	"Rocky" Style – Alternating Right to Left Foot

Group #7

$\frac{\text{Group } \pi r}{r}$	
ENERGI Station 5 Core & Torso, and Balance Disc	
1	V-Sit
2	Plank with Alternating Knee Tuck

Group #8

Jump Rope Intervals	
1	Two Foot (quickly for 20 then slow down for 10 before repeating)
2	"Rocky" Style – Alternating Right to Left Foot

Group #9

	ENERGI Station 2 Pulling & Back, and Angle Bar & Hurdle	
1	Parallel Bar, Horizontal Pull-Up	
2	Hurdle Dips with Knee Up (holding medicine ball)	

Group #10

Jur	np Rope Intervals
1	Two Foot – side-to-side quickly 20 then slow down for 10 then repeat
2	"Rocky" Style – alternating right to left foot

Group #11

ENERGI Station 2 Angle Bar & Hurdle	
1	Travel Push-Up
2	Lateral Walk Over (holding medicine ball)

Group #12

Jump Rope Intervals	
1	Two foot – front to back as quick as you can go for 20 then slow down for 10 then repeat
2	Rocky Style – alternating right to left foot

Cool down & stretch together.

Additional materials needed: lbs. Activity Cards (See pg. 4 for more information) 12 Jump Ropes 2 Medicine Balls